
Samaris™

PEEK ANTERIOR ANATOMICAL
CERVICAL CAGE


SURGICAL TECHNIQUE GUIDE

Scient'x® *Alpha* **Alphatec Spine®**
Solutions for the Aging Spine®

Perfect anatomical design for optimum stability

Optimized fusion

Wide open design to optimize graft area and to enhance bone fusion.

Excellent stability

Inferior and superior **teeth** on surface designed to prevent migration.





Original design to maximize contact with the cortical parts of the vertebral endplates for optimum primary stability.

Anatomical shape / Reinforced security

Biconvexity in the frontal and sagittal planes designed to fit to the intersomatic space

Various footprints and heights designed to adapt to any morphology, restore interbody height and preserve cervical lordosis.

Large cage to support load distribution and optimal attachment on the cortical bone.

Interbody arthrodesis is a secure and efficient technique for the treatment of cervical spine pathologies.

The Samarys™ PEEK anterior anatomical cervical cage fulfills the requirements for this technique.

Its unique design allows for optimum interbody fusion. The anatomical shape and the pure PEEK composition of this implant help to obtain excellent bone fusion. Its broad range of footprints and heights permit to answer the needs of all situations and anatomical variability.

The Samarys™ anterior cervical cage is a surgical implant indicated for cervical arthrodesis via anterior approach, designed to achieve bone fusion from C3 to C7 levels.

It is an interbody cage whose function is to maintain the discal height and the radicular decompression by interbody distraction. Its design facilitates osteogenesis.

The indications are:

- Stabilization after treatment of a herniated cervical intervertebral disc,
- Osteophytosis compressing the nerve roots and/or spinal cord,
- Degenerative intervertebral instability.

Features and Benefits of PEEK-OPTIMA®

PEEK-OPTIMA® by Invibio is a pure PEEK material designed to meet spine, hip and dental implant requirements.

Safety

- Elasticity designed to prevent subsidence: decreased risk of endplate penetration
- Rupture modulus is lower than metal: decreased risk of bone damage
- Allows for repeated autoclave sterilizations

Medical imaging (X-ray, CT, MRI)


- Radiolucent material to allow fusion follow-up
- Embedded gold markers to facilitate implant placement verification

Biocompatibility

- No carbon fibers: decreased risk of inflammatory reaction
- Long history of implantation (since 1999) and successful use in regulatory approved devices (CE marked and FDA cleared devices)

Enhanced bone fusion

- Load sharing effect
- Optimum load repartition
- Bone growth enhancement due to micromovement


Bibliography:

- Pascal-Mousselard H., Daculsi G., Lazennec JY, Catonné Y, Saillant G., Les cages intervertébrales cervicales: analyse critique de la littérature, Maîtrise Orthopédique n°147 - Oct 2005.
- Pascal-Mousselard h., Anterior Cervical Fusion with PEEK Cages: Clinical Results of a Prospective, Comparative, Multicenter and Randomized Study Comparing Iliac Graft and a Macroporous Biphasic Calcium Phosphate, The Spine Journal, 2006, Volume 6, Issue 5, Pages 136S-136S.

References obtained from congress proceedings:

- Ganem F, Bernard P, Adam Y, Hansen F. Cervical Grafts in Degenerative Disease. Poster session presented at 18th annual meeting Cervical Spine Research Society; 2002 June 13-14; Paris, France.
- Bernard P, Ganem F, Ferreira A. Comparison of a Cervical Fusion Cage with Autologous Bone and with Biphasic Ceramic Bone Substitute. Proceeding 2003 June 19-20; Barcelona, Spain.
- Pascal-Mousselard H, Aguado E, Catonné Y, Rouvillain J.L, Daculsi G., Evaluation of an intervertebral cervical cage filled with injectable bone substitute. Proceedings of the SpineWeek International Annual Congress. Processing 2004 May 30 - June 5; Porto, Portugal.


Surgical protocol/patient positioning

The procedure is performed under general anesthesia with endotracheal intubation.

The patient is positioned in the supine position on a standard table. The head is placed in a neutral position or rotated in the opposite direction to the approach. Rotation should be approximately 30° to limit stress on the sterno-cleido-mastoid muscle. The shoulders are lowered and fixed with self-adhesive strips in order to release the inferior cervical spine.


The skin incision targets the level of the lesion to be treated. It is marked on the skin after locating the lesion with an image intensifier. It may be horizontal, along a skin fold if it is well located. It may also be vertical along the sterno-cleido-mastoid for more than two levels.

The platysma muscle is sectioned on the same plane as the skin; then two upper and lower detachments are created to avoid muscle tension.

The aponeurotic planes are sectioned with the finger up to the anterior surface of the spine.

The muscles along the neck are separated and rugined to thoroughly release the operating area and allow proper retraction using an autostatic distractor.

Under scopic control, the correct level is located.

The use of a Caspar distractor might be useful: first the cervical distraction pins are positioned using the screwdriver for distraction pin.

They should be inserted in the middle of the vertebral bodies' over- and under-lying the discal space to be treated and parallel to the endplates.

All the threads of the pins should be engaged in the vertebral bodies. After removing the screwdriver, the cervical distractor can be slid onto the distraction pins until it reaches the base of the pins (closest to the vertebral body).

Then the anterior vertebral ligament is incised down to the lateral unco-vertebral articular processes before distracting the space.

Note: In case of further utilization of the trial/cage holder with stopper, in a shallow intervertebral space, do not forget to position the cervical distraction pins as to let enough space for further holder with stopper introduction.

Discectomy and preparation of the surgical site

The discectomy starts with the rectangular excision of the anterior part of the annulus fibrosus using a scalpel.

Distraction pin Ø 3.5 mm
21FCD35-xx

Screwdriver for distraction pin
21TRV03

Distractor
21DST02

The disc material is resected using curettes and rongeurs. After full discectomy, resection of the posterior osteophytes and the posterior part of the uncus may provide foramen release. For disc herniations, it may be useful to open the posterior longitudinal ligament to ensure all disc material has been removed. It may prove necessary to use a microscope for resecting the posterior disc material.

Preparing the vertebral endplates

The vertebral plates should be prepared and cleaned carefully with a curette. A damaged surface may cause subsidence of the cage in the vertebral body. This procedure is essential to ensure optimum contact between the implant and adjacent vertebral endplates and to avoid secondary migration. The removal of the subchondral parts also improves the bone ingrowth through the graft. In order to help you in this disc removal step, a **cervical rasp** (height: 4,5mm) reproducing the shape of the Samarys cage is also included in the Samarys set.

Option: Several cervical rasps are also available upon request in 6 additional heights (ranging from 5 to 10mm). The rasps should be selected according to the intervertebral height and carefully inserted until resistance is felt.

Note: It is important to use these instruments with caution so as not to excessively rasp the endplates and to insert the instrument too posteriorly.

The following steps describe the surgical technique without the distractor.



Insertion of the trial implants

A set of **trial implants** are used to determine the ideal implant size. Trial implants are color coded according to the footprints, facilitating the identification of the corresponding cage. The heights can be checked by the laser-marking. They range from 4.5mm to 10mm.

The trial is secured onto the **cage holder** by tightening the inner part of the instrument into the central hole of the trial. A horizontal groove improves the stability during impaction. It is then impacted into the intervertebral space. Having no retaining edges, the trial can be easily removed.

Note: To determine the size of the most suitable cage, it is recommended to first insert the smallest trial implant, and sequentially select the next bigger one till the proper size is achieved. It is recommended to use the biggest implant possible considering patient anatomy. Choice of the adequate size is then left to the surgeon's discretion. The height to be restored can be determined according to the adjacent disc heights or the facet parallelism at the index level.

Note: If caspar distraction is used, it is mandatory to release distraction after insertion of the trial.

Trial implant
 21CC2A-15-xxF
 21CC2AC-17-xxF
 21CC2A-17-xxF
 21CC2A-20-xxF

Cage holder
 21ICC2A01

Cervical Rasp
 21RAP01-xx

Graft filling

Bone graft material (autologous, bone substitute...) has to be compacted into the implant. To help this compaction, the cage should be positioned on the cage **socket**.

The upper and lower faces of the cage have wide openings to offer optimum graft volume and contact area. A **graft pusher** matching the internal cavity shape of the cage is used to compact the graft (two different graft pushers are available and have to be chosen according to the footprint of the selected cage).

Cage insertion

After filling, the cage is carefully secured onto the **cage holder** or onto the **cage holder with stopper**.

To hold the cage, the inner part of the holder is fixed into the front hole of the cage (same maneuver as for the trial). Be careful to not too excessively tighten the holder onto the cage as to preserve the thread and PEEK cage integrity). Double check the cage orientation (sagittal curvature at the top). To help, the cage also features a triangle marking (Δ) on its lateral face that allows to identify the top of the cage. The holder has small dimensions that prevent visual obstruction of the implant.

The cage should be carefully impacted with the distractor slightly distracted in order to allow for introduction of the cage while preventing the risk to impact it too posteriorly. When using the cage holder with the stopper, the circular stop at the extremity of the instrument prevents the cage from moving too far posteriorly. The anterior edge of the cage stops at 1.5mm from the anterior wall of the vertebrae. Caution: The cages from 8mm have a height superior to that of the stopper; hence the stop is not functional for cages with a height of 8mm and higher. Although the cage holder features a stopper, it is advisable to insert the cage with caution and vigilance. Introducing the cage under C-arm control can be helpful in this respect.

The holder is removed after implanting the cage by releasing its inner part.

Note: The Samarys cervical cage ensures immediate post-operative stability of the spine but, in some cases, internal fixation, cervical plate,... or external collar brace, can be added at the surgeon's discretion.

Closing the approach wound

The approach wound is closed after rinsing. The haemostasis is checked and a suction drain should be placed on the anterior face of the spine. The platysma muscle must be carefully restored, and the skin closed by intradermal stitches with resorbable thread.

Removing the implant (optional)

If the instrumentation must be removed, before fusion, the cervical approach is used, down to the instrumented area. Bone bridges between the implant and the vertebral bodies must be sectioned. The **cage holder** is placed onto the cage to facilitate implant removal.

Socket
21SOC02


Cage holder with stopper
21ICC2A03

Graft pusher
21COM03
21COM07

Implants

Teeth height: 0.5 mm

PEEK Anterior Anatomical Cervical Cage

Height (mm)	Depth 12 mm Width 15 mm	Depth 13 mm Width 17 mm	Depth 15 mm	
	Ref.	Ref.	Width 17 mm Ref.	Width 20 mm Ref.
4.5	11CC2A15-45S*	11CC2AC17-45S*	11CC2A17-45S*	11CC2A20-45S*
5	11CC2A15-5S	11CC2AC17-5S	11CC2A17-5S	11CC2A20-5S*
6	11CC2A15-6S	11CC2AC17-6S	11CC2A17-6S	11CC2A20-6S*
7	11CC2A15-7S*	11CC2AC17-7S*	11CC2A17-7S	11CC2A20-7S*
8	11CC2A15-8S*	11CC2AC17-8S*	11CC2A17-8S*	11CC2A20-8S*
9	11CC2A15-9S*	11CC2AC17-9S*	11CC2A17-9S*	11CC2A20-9S*
10	11CC2A15-10S*	11CC2AC17-10S*	11CC2A17-10S*	11CC2A20-10S*

The implants are delivered sterile.

Instruments

Trial implants for:

	Color code	Ref.
11CC2A15-xx	Gold	21CC2A15-xxF
11CC2AC17-xx	Silver	21CC2AC17-xxF
11CC2A17-xx	Green	21CC2A17-xxF
11CC2A20-xx	Fuschia	21CC2A20-xxF

NOTE: Trials corresponding to implants size on request, are also available upon request. Delivery time to be confirmed with order.

Cervical distraction pin Ø 3.5 mm
21FCD35-xx

Length (mm)	Ref.
12	21FCD35-12*
14	21FCD35-14
16	21FCD35-16

Cervical distractor 21DST02

Screwdriver for distraction pin 21TRV03

21RAP01-5* 21RAP01-8*
21RAP01-6* 21RAP01-9*
21RAP01-7* 21RAP01-10*
21RAP01-45

Cage holder 21CC2A01

Cage holder with stopper 21CC2A03

Graft pusher 21COM03
21COM07

Socket 21SOC02

*Upon request. Delivery time to be confirmed with order.

! Please refer to the Instructions For Use document enclosed within product packaging prior to use.

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The **SAMARYS anterior cervical cage** is an implant intended for the surgical treatment of the cervical vertebrae.

GENERAL DESCRIPTION:

The **SAMARYS anterior cervical cage** has been designed for stabilisation and arthrodesis between the cervical vertebrae.

The **SAMARYS anterior cervical cage** is convex in shape on its upper face and has openings on the upper and bottom faces.

It comes in several heights to adapt to the variable morphology of the intervertebral spaces.

The **SAMARYS anterior cervical cage** is made of implantable PEEK and includes radiopaque markers made of 99.99% pure gold.

The **SAMARYS anterior cervical cage** must not be used with components coming from other manufacturers.

INDICATIONS:

The **SAMARYS anterior cervical cage** is a surgical implant for cervical arthrodesis by anterior approach, designed to optimize bone fusion from C3 to C7 levels. It is an intersomatic cage whose main function is to maintain the discal height and the radicular decompression by intersomatic distraction. Its design is favouring osteogenesis. The indications are:

- stabilisation after treatment of a herniated cervical intervertebral disc or osteophytosis compressing the nerve roots and/or spinal cord,
- degenerative intervertebral instability.

CONTRA-INDICATIONS:

The contra-indications for the **SAMARYS anterior cervical cage** include:

- local infection or inflammation,
- vertebral osteoporosis,
- malignant vertebral disease,
- allergy or intolerance to PEEK,
- Incompatible patient age and physical condition,
- all cases not included in the indications.

The **SAMARYS anterior cervical cage** has not been designed, intended or sold for uses other than those indicated.

POSSIBLE ADVERSE EFFECTS:

- Infection,
- pain,
- pseudarthrosis,
- adjacent segment disease,
- post-operative migration of the implant ,
- damage to the vertebrae adjacent to the arthrodesis,
- Intolerance to the material.

Note: An additional surgical operation may be needed to correct any adverse event.

Warnings: An entirely satisfactory result is not always obtained at each and every operation. This is particularly true in spinal surgery where many external factors can compromise the results.

OPERATING PRECAUTIONS:

The surgeon must be thoroughly familiar with the **SAMARYS anterior cervical cage**, the method of application, the instruments and the operating technique.

The size of the **SAMARYS anterior cervical cage** must be chosen in relation to the individual clinical case and the desired correction.

It is recommended to use a Caspar type distractor to open the intervertebral space intended for implanting the cage.

Before implanting the **SAMARYS anterior cervical cage**, a careful curettage of vertebral endplates must be conducted to expose healthy tissue but without weakening them to prevent the cage from subsidence.

The **SAMARYS anterior cervical cage** must be filled with autologous or allogenic bone or with a bone substitute to obtain bone fusion.

The correct positioning of the cage with respect to the vertebrae is confirmed by radiographic control during surgery.

After implantation, the **batch number** and the **reference of the implanted SAMARYS anterior cervical cage** must always be recorded in the patient's surgical records.

This product is a single use device. Under no circumstances should it be reused. While the device may appear to be undamaged, it may have small defects or internal stress patterns, as a result of the prior implantation or removal that could lead to fatigue failure. Additionally, please note that the removed implant has not been

designed or validated so as to allow for decontamination of microorganisms. Reuse of this product could lead to cross-infection and/or material degradation as a result of the decontamination process. The company accepts no responsibility for products which have been reused.

Caution: no silicone, formaldehyde or latex based products should be implanted with a SCIENT'X implant.

Patients must be informed of the precautions to be taken in their everyday life to guarantee a maximum implant service life. It is recommended that a regular postoperative follow-up is undertaken to detect early signs of failure of the implants and to consider the action to be taken. Deterioration of the device after bone consolidation cannot be considered to constitute a dysfunction or deterioration in the characteristics of the implant. A suitable rehabilitation program must be designed and implemented.

PACKAGING:

The storage conditions must enable implants, the associated instruments and their respective packages integrity to be maintained.

The good condition of the implants and the functionality of the instruments must be checked before use. Do not expose the SAMARYS implants to radiations or extreme temperature. Should these requirements not be followed, reduced mechanical properties may occur which could lead to implant failure in some cases. Damaged products must never be used and must be returned to SCIENT'X. Do not implant the cage if the primary and/or the secondary package of the sterile conditioning are pierced or deteriorated.

DECONTAMINATION, CLEANING AND STERILIZATION:

Products delivered in sachets are not sterile.

For implants delivered sterile: the implants are sterilized by Gamma radiation at doses of 25 to 40 kGy. The expiry date is 5 years. The expiry date of sterile parts is indicated on the packaging.

For implants and instruments not sterile on delivery: all implants and instruments delivered non-sterile must be **decontaminated, cleaned and sterilized** before and after use. Implants and instruments in sachets must be removed from the original packaging for the following operations:

Recommended method:

- **Decontamination:** Plunge the implants and instruments into a bactericidal and fungicidal solution (didecyl ammonium chloride with proteolytic enzymes diluted to 0.5 % (5g for 1 litre of mild water). Deep length: 30 min. Rinse with demineralized water.
- **Cleaning:** Wash the implants and instruments in a LANCER type machine with suitable cleaning products, rinse and dry. Any product that might damage the equipment is forbidden (such as bleach, formol, etc.). Lubricate all hinge and mechanism on instruments with paraffin oil (type Neodisher IP spray).
- **Sterilization:** It is essential to sterilise the kit by steam using the following conditions:
 - pre-heating for 25' at 110° (1 bar)
 - vacuum 5' (0.8 bar under atmospheric pressure)
 - heating 5' at 120° (1 bar)
 - vacuum 5' (0.8 bar)
 - sterilization 18' at 134° (2 bars)
 - drying 20' return to room temperature

COMPLAINTS:

Any complaints, together with the reference and lot number of the incriminated product, must be sent to SCIENT'X.


ADDITIONAL INFORMATION:

For further information, contact the manufacturer:

SCIENT'X
22, rue Jean Bart
78960 Voisins le Bretonneux - FRANCE
tel.: 33 (0)1 39 30 69 30
fax: 33 (0)1 30 43 82 77
e-mail: qualite@scientx.com

CE0459 2001

Last revision: 06/2010

Only for information.

Subject to modification without notification.

For updated IFU, please contact your local customer service.

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Carlsbad, CA 92008 - U.S.A.